

## TRANSFER BLOOD RESULTS TO LAPTOP



### Steps:

1. Blood machine (**analyzer**) is **Off**
2. Plug **printer cable** into **Laptop** and **Analyzer**.
3. Turn **analyzer On**, allow it to boot to **menu**.
- 3a. Ensure the **analyzer** is **not running** a **rotor** analysis.
4. Run **Abaxis Data Manager** through the **Laptop**.
- 4a. Allow **Abaxis Data Manager** to automatically **query** **analyzers**.
- 4b. If **Abaxis Data Manager** failed to **automatically query**, then click **file**, and **query** **analyzers**.
5. From the **Blood Machine** Home Page **click** the **folder icon**
- 5a. Select results by either **Browsing**, **Searching**, or by **Date Range**
6. **Transmit** results from **Blood Machine** to Laptop.
- 6a. Wait for the transmitting to complete.
7. On the Laptop run the **Abaxis XML Reader**
8. **Choose File** and click **convert**, alternatively you can convert all.

